

Rooke® | VASCULAR MITT

Rooke®
Helping Limbs Heal

RookeProducts.com | 800•535•5865

Instructions for Use:

1. Insert hand so grip material is on the palm-side of the hand



2. Pull the strap for a snug fit



3. Secure the strap



Washing instructions: Wash in cold water 30C/85F, gentle/delicate cycle, line dry, do not bleach, do not dry clean, and do not iron.

- Follow all written instructions included with the product, along with all verbal instructions and advice provided by your physician or care provider.
- **Each Rooke Vascular Mitt Includes:**
One Rooke Vascular Mitt **RV-51**



Rooke | Osborn Medical Corporation

9800 E Easter Ave, Suite 130 Centennial, CO 80112

Call: USADirect: 303•223•1800 | Toll Free: 1•800•535•5865 **Fax:** 507•932•5044

For more information about the Rooke Vascular Mitt please visit: www.RookeProducts.com