

Care Instructions

- Washing instructions: Wash in cold water 30C/85F, gentle/delicate cycle
- Drying instructions: Line dry or tumble dry on low temperature
- DO NOT: bleach, dry clean, or iron.

Rooke | Osborn Medical Corporation

9800 E Easter Ave, Suite 130 Centennial, CO 80112

Call: USADirect: 303•223•1800 | Toll Free: 1•800•535•5865 **Fax:** 507•932•5044

For more information and detailed videos of the Rooke Vascular Boot please visit: rookeproducts.com